




ECTON (Sine Wave)



Features :

1. Digital Display of voltage.
2. Voltage selector from 80 to 180 v at 10 v step.
3. Four models of operation.
4. Manual, Auto, Ectonus & CST Mode.
5. For Auto-mode timer from 0 to 1.6 sec.
6. Ectonus through liner Variant.
7. CST pulses from 0 to 40 v to reduce amnesia & confusion after ECT & to revive memory. It is a non convulsive treatment & is beneficial in NEUROTIC cases too.
8. Highly efficient spring loaded Head band to avoid skin burning.
9. Treat switch provided on Headband too.
10. Ectonus gives the smooth Onset of the convulsive stage, controls the colonic stage & avoids the need of anesthetics in all normal cases.
11. Velcro tape electrodes for CST.
12. Transformer fully isolated.
13. Odourless mouth gags.
14. Mains operated at 230 v 50 Hz
15. Compatible to Foot switch & unilateral electrodes.

RECORDERS & MEDICARE SYSTEMS (P) LTD.

181/5, Industrial Area, Phase-1, Chandigarh 160 002, INDIA.
Ph. : 91-172-2658701-705 Fax : 91-172-2653415, 5073626
e-mail : rmshelpdesk@rmsindia.com
website : www.rmsindia.com

An ISO 9001 : 2000 company